

MT-ND1 Antibody



PACO13533

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MT-ND1

Uniprot

P03886

Synonyms:

NADH dehydrogenase, subunit 1 (complex I); ND1; MTND1 ;

Immunogen:

Human ND1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A