

MT-ND2 Antibody



PACO13534

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MT-ND2

Uniprot

P03891

Synonyms:

NADH dehydrogenase, subunit 2 (complex I); ND2; MTND2 ;

Immunogen:

Human ND2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A