

PACO13538

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GNAS

Uniprot

O95467

Synonyms:

GNAS complex locus; GNAS; AHO; C20orf45; GNAS1; GPSA; GSA; GSP; MGC33735; NESP; PHP1A; PHP1B; POH; DJ309F20.1.1; DJ806M20.3.3 ;

Immunogen:

Human NESP55, GNAS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A