

# ILF2 Antibody



PACO13539

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ILF2

**Uniprot**

Q12905

**Synonyms:**

interleukin enhancer binding factor 2, 45kDa; ILF2; MGC8391; NF45; PRO3063 ;

**Immunogen:**

Human NF45.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A