

ILF3 Antibody



PACO13540

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ILF3

Uniprot

Q12906

Synonyms:

interleukin enhancer binding factor 3, 90kDa; ILF3; CBTF; DRBF; DRBP76; MMP4; MPHOSPH4; MPP4; NF-AT-90; NF110; NF90; NFAR; NFAR-1; NFAR2; TCP110; TCP80 ;

Immunogen:

Human NF90, ILF3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A