

NPRL2 Antibody



PACO13543

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NPRL2

Uniprot

Q8WTW4

Synonyms:

44kd;

Immunogen:

Human NPR2L, TUSC4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A