

# IST1 Antibody



PACO13548

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IST1

**Uniprot**

P53990

**Synonyms:**

KIAA0174; KIAA0174; IST1; MGC117220 ; OLC1

**Immunogen:**

Human OLC1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A