

CDKN2A Antibody



PACO13551

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CDKN2A

Uniprot

P42771

Synonyms:

cyclin-dependent kinase inhibitor 2A (melanoma, p16, inhibits CDK4); CDKN2A; ARF; CDK4I; CDKN2; CMM2; INK4; INK4a; MLM; MTS1; TP16; p14; p14ARF; p16; p16INK4; p16INK4a; p19 ; ARF, CDKN2A

Immunogen:

Human P16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A