

DDX6 Antibody



PACO13552

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DDX6

Uniprot

P26196

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 6; DDX6; FLJ36338; HLR2; P54; RCK ;

Immunogen:

Human P54.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A