

PADI2 Antibody



PACO13557

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PADI2

Uniprot

Q9Y2J8

Synonyms:

peptidyl arginine deiminase, type II; PADI2; KIAA0994; PAD-H19; PAD2; PDI2 ;

Immunogen:

Human PADI2, PAD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A