

PARK7 Antibody



PACO13558

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PARK7

Uniprot

Q99497

Synonyms:

Parkinson disease (autosomal recessive, early onset) 7; PARK7; DJ-1; DJ1; FLJ27376; FLJ34360; FLJ92274 ;

Immunogen:

Human PARK7, DJ-1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A