

POLR2B Antibody



PACO13567

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

POLR2B

Uniprot

P30876

Synonyms:

polymerase (RNA) II (DNA directed) polypeptide B, 140kDa; POLR2B; POL2RB; RPB2; hRPB140; hsRPB2 ;

Immunogen:

Human POLR2B-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A