

POTED Antibody



PACO13568

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POTED

Uniprot

Q86YR6

Synonyms:

POTE ankyrin domain family, member D; POTED; A26B3; ANKRD21; POTE; POTE-21; POTE21 ;

Immunogen:

Human POTED-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A