

# POU3F3 Antibody



PACO13570

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

POU3F3

**Uniprot**

P20264

**Synonyms:**

POU class 3 homeobox 3; POU3F3; BRN1; OTF8 ;

**Immunogen:**

Human POU3F3-Specific.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A