

PRDX5 Antibody



PACO13576

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PRDX5

Uniprot

P30044

Synonyms:

peroxiredoxin 5; PRDX5; ACR1; AOEB166; B166; MGC117264; MGC142283; MGC142285; PLP; PMP20; PRDX6; PRXV; SBBI10 ;

Immunogen:

Human Prx5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A