

PSEN1 Antibody



PACO13577

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PSEN1

Uniprot

P49768

Synonyms:

presenilin 1; PSEN1; AD3; FAD; PS1; S182 ;

Immunogen:

Human PSEN1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A