

REXO4 Antibody



PACO13589

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

REXO4

Uniprot

Q9GZR2

Synonyms:

REX4, RNA exonuclease 4 homolog (S. cerevisiae); REX4; REX4; XPMC2; XPMC2H ;

Immunogen:

Human REXO4, REX4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A