

RIPK1 Antibody



PACO13594

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RIPK1

Uniprot

Q13546

Synonyms:

receptor (TNFRSF)-interacting serine-threonine kinase 1; RIPK1; FLJ39204; RIP; RIP1 ;

Immunogen:

Human RIPK1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A