

# ROBO4 Antibody



PACO13596

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ROBO4

**Uniprot**

Q8WZ75

**Synonyms:**

roundabout homolog 4, magic roundabout (Drosophila); ROBO4; ECSM4; FLJ20798; MGC133352; MGC133353; MRB ;

**Immunogen:**

Human ROBO4-Specific.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A