

RPS4X Antibody



PACO13598

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RPS4X

Uniprot

P62701

Synonyms:

ribosomal protein S4, X-linked; RPS4X; CCG2; DXS306; FLJ40595; SCAR; SCR10 ;

Immunogen:

Human SCAR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A