

OXCT1 Antibody



PACO13600

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

OXCT1

Uniprot

P55809

Synonyms:

3-oxoacid, CoA transferase 1; OXCT1; OXCT; SCOT ;

Immunogen:

Human SCOT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A