

SDC1 Antibody



PACO13602

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SDC1

Uniprot

P18827

Synonyms:

syndecan 1; SDC1; CD138; SDC; SYND1; syndecan ;

Immunogen:

Human SDC1, CD138.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A