

SENP5 Antibody



PACO13604

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SENP5

Uniprot

Q96HI0

Synonyms:

SUMO1/sentrin specific peptidase 5; SENP5; DKFZp564O1016; FLJ42398; MGC27076 ;

Immunogen:

Human SENP5-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A