

SEPT2 Antibody



PACO13607

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

37500

Uniprot

Q15019

Synonyms:

septin 2; SEPT2; DIFF6; KIAA0158; NEDD5; Pnutl3; hNedd5 ; septin2

Immunogen:

Human SEPT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A