

SEPT9 Antibody



PACO13609

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

40057

Uniprot

Q9UHD8

Synonyms:

septin 9; SEPT9; AF17q25; FLJ75490; KIAA0991; MSF; MSF1; NAPB; PNUTL4; SINT1; SeptD1 ;

Immunogen:

Human SEPT9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A