

ATP2A2 Antibody



PACO13610

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATP2A2

Uniprot

P16615

Synonyms:

ATPase, Ca++ transporting, cardiac muscle, slow twitch 2; ATP2A2; ATP2B; DAR; DD; DKFZp686P0211; FLJ20293; FLJ38063; MGC45367; SERCA2 ;

Immunogen:

Human SERCA2, ATP2A2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A