

SERINC2 Antibody



PACO13611

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SERINC2

Uniprot

Q96SA4

Synonyms:

serine incorporator 2; SERINC2; FKSG84; MGC90340; PRO0899; TDE2; TDE2L ;

Immunogen:

Human SERINC2-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A