

SERPINA3 Antibody



PACO13612

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SERPINA3

Uniprot

P01011

Synonyms:

serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 3; SERPINA3; AACT; ACT; GIG24; GIG25; MGC88254 ; A1-Antichymotrypsin

Immunogen:

Human SERPINA3, AACT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A