

# SKAP1 Antibody



PACO13614

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SKAP1

**Uniprot**

Q86WW1

**Synonyms:**

src kinase associated phosphoprotein 1; SKAP1; SCAP1; SKAP55 ;

**Immunogen:**

Human SKAP55, SKAP1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A