

SLC12A5 Antibody



PACO13615

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC12A5

Uniprot

Q9H2X9

Synonyms:

solute carrier family 12 (potassium-chloride transporter), member 5; SLC12A5; KCC2; KIAA1176 ;

Immunogen:

Human SLC12A5-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A