

# SLC2A1 Antibody



PACO13616

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC2A1

**Uniprot**

P11166

**Synonyms:**

solute carrier family 2 (facilitated glucose transporter), member 1; SLC2A1; DYT17; DYT18; GLUT; GLUT1; MGC141895; MGC141896; PED ;

**Immunogen:**

Human SLC2A1, GLUT1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A