

SMTNL1 Antibody



PACO13621

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SMTNL1

Uniprot

A8MU46

Synonyms:

smoothelin-like 1; Smtnl1;1110030K22Rik; CHASM ;

Immunogen:

Human Smtnl1-phospho.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A