

SMYD1 Antibody



PACO13622

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SMYD1

Uniprot

Q8NB12

Synonyms:

SET and MYND domain containing 1; SMYD1; BOP; ZMYND18; ZMYND22 ;

Immunogen:

Human SMYD1-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A