

ST8SIA2 Antibody



PACO13628

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ST8SIA2

Uniprot

Q92186

Synonyms:

ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 2; ST8SIA2; HsT19690; MGC116854; MGC116857; SIAT8B; ST8SIA-II; STX ;

Immunogen:

Human ST8SIA2-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A