

STEAP2 Antibody



PACO13629

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

STEAP2

Uniprot

Q8NFT2

Synonyms:

six transmembrane epithelial antigen of the prostate 2; STEAP2; IPCA1; PCANAP1; PUMPCn; STAMP1; STMP ;

Immunogen:

Human STEAP2-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A