

BIRC5 Antibody



PACO13630

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BIRC5

Uniprot

O15392

Synonyms:

baculoviral IAP repeat-containing 5; BIRC5; API4; EPR-1 ; ; SURVIVIN; BIRC5

Immunogen:

Human SURVIVIN.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A