

TFAP2A Antibody



PACO13633

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TFAP2A

Uniprot

P05549

Synonyms:

transcription factor AP-2 alpha (activating enhancer binding protein 2 alpha); TFAP2A; AP-2; AP-2alpha; AP2TF; BOFS; FLJ51761; TFAP2 ; AP2

Immunogen:

Human TFAP2A, AP-2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A