

TGFBR3 Antibody



PACO13635

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TGFBR3

Uniprot

Q03167

Synonyms:

transforming growth factor, beta receptor III; TGFBR3; BGCAN; betaglycan ; TGF RIII

Immunogen:

Human TGFBR3-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A