

# PYCARD Antibody



PACO13641

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PYCARD

**Uniprot**

Q9ULZ3

**Synonyms:**

PYD and CARD domain containing; PYCARD; ASC; CARD5; MGC10332; TMS; TMS-1; TMS1 ;

**Immunogen:**

Human TMS1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A