

TNFRSF1A Antibody



PACO13642

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TNFRSF1A

Uniprot

P19438

Synonyms:

tumor necrosis factor receptor superfamily, member 1A; TNFRSF1A; CD120a; FPF; MGC19588; TBP1; TNF-R; TNF-R-I; TNF-R55; TNFAR; TNFR1; TNFR55; TNFR60; p55; p55-R; p60 ;

Immunogen:

Human TNFR1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A