

TARBP2 Antibody



PACO13649

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TARBP2

Uniprot

Q15633

Synonyms:

TAR (HIV-1) RNA binding protein 2; TARBP2; LOQS; TRBP; TRBP1; TRBP2 ;

Immunogen:

Human TRBP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A