

TSPYL2 Antibody



PACO13652

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TSPYL2

Uniprot

Q9H2G4

Synonyms:

TSPY-like 2; TSPYL2; CDA1; CINAP; CTCL; DENTT; HRIHFB2216; SE20-4 ;

Immunogen:

Human TSPYL2, TSPX, NP79.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A