

U2AF1 Antibody



PACO13655

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

U2AF1

Uniprot

Q01081

Synonyms:

U2 small nuclear RNA auxiliary factor 1; U2AF1; DKFZp313J1712; FP793; RN; RNU2AF1; U2AF35; U2AFBP ;

Immunogen:

Human U2AF35.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A