

# URB1 Antibody



PACO13663

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

URB1

**Uniprot**

O60287

**Synonyms:**

URB1 ribosome biogenesis 1 homolog (S. cerevisiae); URB1; C21orf108; KIAA0539; NPA1 ;

**Immunogen:**

Human URB1-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A