

PACO13665

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GC

Uniprot

P02774

Synonyms:

group-specific component (vitamin D binding protein); GC; DBP; DBP/GC; VDBG; VDBP ;

Immunogen:

Human VDBP, GC.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A