

SLC18A2 Antibody



PACO13666

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC18A2

Uniprot

Q05940

Synonyms:

solute carrier family 18 (vesicular monoamine), member 2; SLC18A2; MGC120477; MGC120478; MGC26538; SVAT; SVMT; VAT2; VMAT2 ;

Immunogen:

Human VMAT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A