

ZIC4 Antibody



PACO13668

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZIC4

Uniprot

Q8N9L1

Synonyms:

Zic family member 4; ZIC4; FLJ42609; FLJ45833 ;

Immunogen:

Human ZIC4-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A