

AQP4 Antibody



PACO13672

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AQP4

Uniprot

P55087

Synonyms:

HMIWC2; MGC22454; MIWC; aquaporin 4

Immunogen:

Human AQP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A