

EDN3 Antibody



PACO13681

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EDN3

Uniprot

P14138

Synonyms:

endothelin 3; EDN3; ET3; MGC15067; MGC61498 ;

Immunogen:

Human EDN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A