

# FZD10 Antibody



PACO13685

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FZD10

**Uniprot**

Q9ULW2

**Synonyms:**

frizzled homolog 10 (Drosophila); FZD10; CD350; FZ-10; FzE7; hFz10 ;

**Immunogen:**

Human FZD10.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A